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Instruction for use

A SOLID-PHASE ENZYME IMMUNOASSAY FOR THE QUANTITATIVE DETERMINATION OF HCG IN HUMAN SERUM OR PLASMA

1. INTENDED USE

A solid-phase enzyme immunoassay for the quantitative determination of hCG in blood serum or plasma.

This kit is designed for measurement of hCG in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Human chorionic gonadotropin (HCG) is a glycoprotein hormone secreted by trophoblastic cells of placenta. A molecule of HCG consists of two noncovalently bound subunits: alpha- and beta-HCG. Beta-subunit is specific for HCG hormone.

Determination of HCG is widely used for early diagnosis of pregnancy. Multiple pregnancy results in correspondent elevation of serum HCG; while ectopic pregnancy and placental insufficiency cause decreased serum HCG levels. Determination of HCG in serum during second trimester is used for pregnancy monitoring, especially in screening for Down syndrome, along with other laboratory tests (AFP and Estriol).

Serum HCG is also a laboratory marker of trophoblastic tumours – chorionepitheliomas, some seminomas and teratomas. Serial determination of serum HCG can be used for therapy monitoring in these cancers.

The present test system uses beta (β)-chain specific monoclonal antibody as the capture, and alpha (α)-chain specific monoclonal antibody as the tracer; therefore only the whole intact HCG molecule is detected.

3. PRINCIPLE OF THE TEST

This test is based on two-site sandwich enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific murine monoclonal to β chain of hCG-antibodies. Antigen from the specimen is captured by the antibodies coated onto the microwell surface. Second antibodies – murine monoclonal to α chain of hCG, labelled with peroxidase enzyme, are then added into the microwells. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is directly related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0% H_2SO_4 . It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

Symbol	Description	Qty	Units	Colour code	Stability of opened/diluted components
1	hCG EIA strips, 8x12 wells	1	pcs		until exp.date
2	Calibrator set, 0.8 ml each. The set contains 6 calibrators: 0, 15, 60, 125, 250, 500 IU/l	6	pcs	blue (CI - colourless)	2 months
3	Control serum (0.8 ml)	1	pcs	colourless	2 months
4	Conjugate, 14 ml	1	pcs	red	until exp.date
5	EIA buffer 100 ml	1	pcs	blue	until exp.date
6	Substrate solution, 14 ml	1	pcs	colourless	until exp.date
7	Washing solution concentrate 26X, 22 ml	1	pcs	colourless	Concentrate - until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT
8	Stop solution, 14 ml	1	pcs	colourless	until exp.date
9	Plate sealing tape	2	pcs		N/A
10	Instruction hCG EIA	1	pcs		N/A
11	QC data sheet hCG EIA	1	pcs		N/A

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for +37 °C ±0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells **TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED)** to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note: It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Alternative units: 1 IU/l = 0.2 ng/ml

7.5. Highdose hook effect. A high-dose hook effect is possible in kits of this nature. To avoid falsely low results, samples obtained from pregnant women should be diluted accordingly before assaying:

Week 1-2 of gestation – no dilution

Week 3-4 – 1:20

Week 4-14 – 1:400

Week 14-21 – 1:400 and 1:1000 (both dilutions should be assayed)

Week 22, 3rd trimester – 1:400.

Dilution examples:

1:20 – 25 µl of sample + 475 µl of EIA buffer

1:400 – 25 µl of 1:20 dilution + 475 µl of EIA buffer

1:1000 – 200 µl of 1:400 dilution + 300 µl of EIA buffer.

If pregnancy status is not clear but suspected, the sample should be assayed both undiluted and diluted 1:20.

7.6. Assay procedure

1	Put the desired number of microstrips into the frame; allocate 14 wells for the calibrators CAL 1–6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS.
2	If suggested analyte concentration in the sample exceeds the highest calibrator, additionally dilute this sample accordingly, using DIL SPE (EIA buffer) (see p 7.5). Use of other buffers or reagents for sample dilution may lead to incorrect measurement.
3	Dispense 100 µl of CONJ HRP into the wells.
4	Pipet 50 µl of calibrators CAL 1–6, control samples CONTROL and unknown samples into the wells. Cover the wells by plate adhesive tape (included into the kit).
5	Incubate 60 minutes at 37 °C.
6	Prepare washing solution by 26X dilution of washing solution concentrate (BUF WASH 26X) by distilled water. Wash the strips 5 times.
7	Dispense 100 µl of SUBS TMB into the wells
8	Incubate 10–20 minutes at +18...+25 °C
9	Dispense 100 µl of STOP into the wells.
10	Measure OD (optical density) at 450 nm.
11	Set photometer blank on first calibrator
12	Apply point-by-point method for data reduction.

7.7. Handling notes

Calibrators and control sample(s) – only one freezing/thawing cycle is allowed

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus hCG concentration.

9.3. Determine the corresponding concentration of hCG in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

Calibrators	Value	Absorbance Units (450 nm)
CAL 1	0 IU/l	0.04
CAL 2	15 IU/l	0.15
CAL 3	60 IU/l	0.41
CAL 4	125 IU/l	0.81
CAL 5	250 IU/l	1.43
CAL 6	500 IU/l	2.36

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for hCG. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

Sex, age	Units, IU/l		Units alternative, ng/ml	
	Lower limit	Upper limit	Lower limit	Upper limit
Males	-	15	-	3.0
Females	-	15	-	3.0
Pregnancy week:				
1	-	50	-	10.0
2	20	500	4.0	100
3	500	5000	100	1000
4	3000	19000	600	3800
5-8	14000	169000	2800	33800
9-13	16000	180000	3200	36000
22	4500	70000	900	14000
23	3000	69500	600	13900
3rd trimester	2400	50000	480	10000

Pregnancy, week	Median, kIU/l	SKO
14	55	0.83
15	42.8	0.78
16	38.4	0.74
17	32.8	0.75
18	25.6	0.69
19	23	0.56
20	21	0.52
21	18	0.51

11. PERFORMANCE CHARACTERISTIC

11.1. Analytical specificity / Cross reactivity

Analyte	Cross-reactivity, % wt/wt
LH	<0.1
FSH	<0.1
TSH	<0.1

11.2. Analytical sensitivity

Sensitivity of the assay was assessed as being 1.25 IU/l.

11.3. Linearity

Linearity was checked by assaying dilution series of 5 samples with different hCG concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery

Recovery was estimated by assaying 5 mixed samples with known hCG concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Krieg, A.F. Pregnancy Tests and Evaluation of Placental Function. Clinical Diagnosis and Management by Laboratory Methods, 16th ed. Henry, J.B. editor, W.B. Saunders Co. Philadelphia. app. 680, 1979.
2. Brody, S., and Carlstrom, G., Immunoassay of Human Chorionic Gonadotropin in Normal and Pathologic Pregnancy. J. Clin. Endocrinol. Metab., 22: 564, 1962.
3. Braunstein, G.D., Rasor, J., Adleer, D., et al: Serum Human Chorionic Gonadotropin Levels Throughout Normal Pregnancy. Am. J. Obstet. Gynecol., 126: 678, 1976.
4. Skinner, M.S., Seckinger, D., Evaluation of Beta-Subunit Chorionic Gonadotropin as an Aid in Diagnosis of Trophoblastic Disease. Ann. Clin. Lab. Sci., 9 (4): 347-52, 1979.
5. Braunstein, G.D., Vaitukaitis, J.L., Carbone, P.P., and Ross, G.T., Ectopic Production of Human Chorionic Gonadotropin by Neoplasms. Ann. Intern. Med., 78: 39-45, 1973.