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*Instruction for use***A SOLID-PHASE ENZYME IMMUNOASSAY
FOR THE QUANTITATIVE DETERMINATION
OF FREE ESTRIOL IN HUMAN SERUM OR PLASMA****1. INTENDED USE**

A solid-phase enzyme immunoassay for the quantitative determination of free estriol in blood serum or plasma.

This kit is designed for measurement of free estriol in blood serum or plasma. For possibility of use with other sample types, please, refer to Application Notes (on request). The kit contains reagents sufficient for 96 determinations and allows to analyze 41 unknown samples in duplicates.

2. SUMMARY AND EXPLANATION

Free (non-conjugated) Estriol (E3) is a steroid hormone with a MM of 288 Da. Production of Estriol in non-pregnant women and men is negligible with no sex-related differences. During pregnancy, Estriol is produced in placenta by a multistep metabolic pathway, and its serum level is continuously growing up to delivery. Due to wide normal ranges of Free E3 at different pregnancy terms, serial determinations to monitor Estriol level are recommended.

Estriol biological function is to maintain blood flow in maternal blood vessels and to provide for differentiation of mammary ducts.

During 2nd trimester, E3 determination in conjunction with other tests (AFP, beta-HCG) is used to estimate risk of Down syndrome in fetus: a tendency of serum E3 level to decline is indicative of Down syndrome. Besides, a decline of E3 serum level during the 2nd trimester may indicate dangerous fetoplacental insufficiency (an increased risk of miscarriage, intra-uterine infections, etc.).

Increased E3 levels are seen in multifetal pregnancy and in case of a fetus bigger than normal.

3. PRINCIPLE OF THE TEST

This test is based on competition enzyme immunoassay principle. Tested specimen is placed into the microwells coated by specific rabbit polyclonal to estriol-antibodies simultaneously with conjugated free Estriol-peroxidase. Free Estriol from the specimen competes with the conjugated free Estriol for coating antibodies. After washing procedure, the remaining enzymatic activity bound to the microwell surface is detected and quantified by addition of chromogen-substrate mixture, stop solution and photometry at 450 nm. Optical density in the microwell is inversely related to the quantity of the measured analyte in the specimen.

4. WARNINGS AND PRECAUTIONS

4.1. For professional use only.

4.2. This kit is intended for in vitro diagnostic use only.

4.3. INFECTION HAZARD: There is no available test methods that can absolutely assure that Hepatitis B and C viruses, HIV-1/2, or other infectious agents are not present in the reagents of this kit. All human products, including patient samples, should be considered potentially infectious. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guidelines or regulations.

4.4. Avoid contact with stop solution containing 5.0% H₂SO₄. It may cause skin irritation and burns.

4.5. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents may give false results.

4.6. Do not use the kit beyond the expiration date.

4.7. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microplate readers.

4.8. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.

4.9. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guidelines or regulations.

4.10. Do not mix reagents from different lots.

4.11. Replace caps on reagents immediately. Do not swap caps.

4.12. Do not pipette reagents by mouth.

4.13. Specimens must not contain any AZIDE compounds – they inhibit activity of peroxidase.

4.14. Material Safety Data Sheet for this product is available upon request directly from XEMA Co., Ltd.

4.15. The Material Safety Data Sheet fit the requirements of EU Guideline 91/155 EC.

5.1. Contents of the Kit

5. KIT COMPONENTS

| Symbol | Description | Qty | Units | Colour code | Stability of opened/ diluted components |
|----------------|--------------------------------------------------------------------------------------------|-----|-------|------------------------|-----------------------------------------------------------------------------------------------|
| 1 SORB MTP | free Estriol EIA strips, 8x12 wells | 1 | pcs | | until exp. date |
| 2 CAL 1-6 | Calibrator set, 0.8 ml each. The set contains 6 calibrators: 0; 0.5; 1.5; 5; 15; 50 nmol/l | 6 | pcs | blue (C1 – colourless) | 2 months |
| 3 CONTROL | Control serum (0.8 ml) | 1 | pcs | colourless | 2 months |
| 4 CONJ HRP | Conjugate, 14 ml | 1 | pcs | green | until exp. date |
| 5 SUBS TMB | Substrate solution, 14 ml | 1 | pcs | colourless | until exp. date |
| 6 BUF WASH 26X | Washing solution concentrate 26X, 22 ml | 1 | pcs | colourless | Concentrate - until exp.date Diluted washing solution - 45 days at 2-8 °C or 15 days at RT |
| 7 STOP | Stop solution, 14 ml | 1 | pcs | colourless | until exp. date |
| 8 N003 | Plate sealing tape | 2 | pcs | | N/A |
| 9 K2181 | Instruction free Estriol EIA | 1 | pcs | | N/A |
| 10 K218Q | QC data sheet free Estriol EIA | 1 | pcs | | N/A |

5.2. Equipment and material required but not provided

- Distilled or deionized water;
- Automatic or semiautomatic multichannel micropipettes, 100–250 µl, is useful but not essential;
- Calibrated micropipettes with variable volume, range volume 25–250 µl;
- Dry thermostat for 37 °C ±0.1 °C;
- Calibrated microplate photometer with 450 nm wavelength and OD measuring range 0–3.0

5.3. Storage and stability of the Kit

Store the whole kit at +2...+8 °C upon receipt until the expiration date.

After opening the pouch keep unused microtiter wells TIGHTLY SEALED BY ADHESIVE TAPE (INCLUDED) to minimize exposure to moisture.

6. SPECIMEN COLLECTION AND STORAGE

This kit is intended for use with serum or plasma (ACD- or heparinized). Grossly hemolytic, lipemic, or turbid samples should be avoided.

Specimens may be stored for up to 48 hours at +2...+8 °C before testing. For a longer storage, the specimens should be frozen at -20 °C or lower. Repeated freezing/thawing should be avoided.

7. TEST PROCEDURE**7.1. Reagent Preparation**

- All reagents (including unsealed microstrips) should be allowed to reach room temperature (+18...+25 °C) before use.
- All reagents should be mixed by gentle inversion or vortexing prior to use. Avoid foam formation.
- It is recommended to spin down shortly the tubes with calibrators on low speed centrifuge.
- Prepare washing solution from the concentrate BUF WASH 26X by 26 dilutions in distilled water.

7.2. Procedural Note:

It is recommended that pipetting of all calibrators and samples should be completed within 3 minutes.

7.3. Assay flowchart

See the example of calibration graphic in Quality Control data sheet.

7.4. Alternative units:

1 nmol/l = 0.288 ng/ml

7.5. Assay procedure

| | |
|----|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1 | Put the desired number of microstrips into the frame; allocate 12 wells for the calibrators CAL 1–6 and control samples CONTROL and two wells for each unknown sample. DO NOT REMOVE ADHESIVE SEALING TAPE FROM UNUSED STRIPS. |
| 2 | Pipet 25 µl of calibrators CAL 1–6, control samples CONTROL and unknown samples into the wells. |
| 3 | Dispense 100 µl of CONJ HRP into the wells. Cover the wells by plate adhesive tape (included into the kit). |
| 4 | Incubate 120 minutes at 37 °C. |
| 5 | Prepare washing solution by 26X dilution of washing solution concentrate (BUF WASH 26X) by distilled water. Wash the strips 5 times. |
| 6 | Dispense 100 µl of SUBS TMB into the wells. |
| 7 | Incubate 10-20 minutes at +18...+25 °C. |
| 8 | Dispense 100 µl of STOP into the wells. |
| 9 | Measure OD (optical density) at 450 nm. |
| 10 | Set photometer blank on air. |
| 11 | Apply lin-log method for data reduction. |

8. QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results.

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable federal, state, and local standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications.

9. CALCULATION OF RESULTS


9.1. Calculate the mean absorbance values (OD450) for each pair of calibrators and samples.

9.2. Plot a calibration curve on graph paper: OD versus free estriol concentration.

9.3. Determine the corresponding concentration of free estriol in unknown samples from the calibration curve. Manual or computerized data reduction is applicable on this stage. Point-by-point or linear data reduction is recommended due to non-linear shape of curve.

9.4. Below is presented a typical example of a standard curve with the XEMA Co. Not for calculations!

| Calibrators | Value | Absorbance Units (450 nm) x 1000 |
|-------------|------------|----------------------------------|
| CAL 1 | 0 nmol/l | 1606 |
| CAL 2 | 0.5 nmol/l | 1268 |
| CAL 3 | 1.5 nmol/l | 1078 |
| CAL 4 | 5 nmol/l | 812 |
| CAL 5 | 15 nmol/l | 569 |
| CAL 6 | 50 nmol/l | 377 |

10. EXPECTED VALUES

Therapeutical consequences should not be based on results of IVD methods alone – all available clinical and laboratory findings should be used by a physician to elaborate therapeutically measures. Each laboratory should establish its own normal range for free Estriol. Based on data obtained by XEMA, the following normal range is recommended (see below). NOTE: the patients that have received murine monoclonal antibodies for radioimaging or immunotherapy develop high titered anti-mouse antibodies (HAMA). The presence of these antibodies may cause false results in the present assay. Sera from HAMA positive patients should be treated with depleting adsorbents before assaying.

| Sex, age | Units nmol/l | | Units alternative ng/ml | |
|------------|-----------------|-------------|----------------------------|-------------|
| | Lower limit | Upper limit | Lower limit | Upper limit |
| Males | - | 1.0 | - | 0.29 |
| Females | - | 1.5 | - | 0.43 |
| Pregnancy: | | | | |
| 22-23 | 5.0 | 20.0 | 1.4 | 5.8 |
| 24-25 | 5.5 | 24.5 | 1.6 | 7.1 |
| 26-27 | 6.3 | 27.9 | 1.8 | 8.0 |
| 28-29 | 6.4 | 31.0 | 1.8 | 8.9 |
| 30-32 | 7.3 | 32.9 | 2.1 | 9.0 |
| 33-34 | 9.9 | 40.3 | 2.9 | 11.6 |
| 35-36 | 12.2 | >50 | 3.5 | >14.4 |
| 37-38 | 15.5 | >50 | 4.5 | >14.4 |
| 39-40 | 16.5 | >50 | 4.8 | >14.4 |
| 40-42 | 17.5 | >50 | 5.0 | >14.4 |

| Pregnancy, week | Median, nmol/l | SKO |
|-----------------|----------------|------|
| 14 | 3.8 | 0.51 |
| 15 | 5.3 | 0.76 |
| 16 | 6.2 | 0.85 |
| 17 | 7.3 | 0.94 |
| 18 | 9.9 | 1.05 |
| 19 | 10.7 | 1.10 |
| 20 | 11.9 | 1.08 |
| 21 | 13.2 | 1.08 |

11. PERFORMANCE CHARACTERISTICS

11.1. Analytical specificity / Cross reactivity

| Analyte | Cross-reactivity, % wt/wt |
|----------------|---------------------------|
| Estradiol | 0.8 |
| Estron | 0.3 |
| 16-epiestriol | 1.9 |
| Corticosterone | <0.1 |

11.2. Analytical sensitivity.

Sensitivity of the assay was assessed as being 0.12 nmol/l.

11.3. Linearity.

Linearity was checked by assaying dilution series of 5 samples with different free estriol concentrations. Linearity percentages obtained ranged within 90 to 110%.

11.4. Recovery. Recovery was estimated by assaying 5 mixed samples with known free estriol concentrations. The recovery percentages ranged from 90 to 110%.

12. LITERATURE

1. Reynold C. Merrill – Estriol: A Review. *Physiol Rev*, Jul 1958; 38: 463 – 480.
2. Karlla K. Welch and Fergal D. Malone – Advances in Prenatal Screening: Maternal Serum Screening for Down Syndrome. *NeoReviews*, Oct 2002; 3: 209.
3. Dan Tulchinsky – Placental Secretion of Unconjugated Estrone, Estradiol and Estriol into the Maternal and the Fetal Circulation. *J. Clin. Endocrinol. Metab.*, Jun 1973; 36: 1079 – 1087.